

Midmark 627 + Ritter 224 Barrier-Free® Examination Chairs

Accessibility in healthcare goes beyond the ability of patients to access the care they need. The physical accessibility of a facility and the equipment within it plays a significant role. The proper exam chair is important in providing a thorough examination. Barrier-Free examination chairs by Midmark are designed to make accessing the exam chair easier, even for those with mobility concerns or other physical limitations.

1. BP POSITIONING

The Midmark® 627 with Patient Support Rails+ helps achieve proper patient positioning for blood pressure capture—arm supported at heart height, back supported and feet flat on the floor.

2. POWERED BACK + HEIGHT*

Powered movement assists in positioning the patient for more comfortable examinations.

3. CLEAN ASSIST ROLLER SYSTEM*

The retractable roller base allows you to easily and safely move the chair for cleaning.

4. STANDARD STIRRUPS

The stirrups are integrated into the exam chair for lower-body examinations.

5. ACTIVE SENSING TECHNOLOGY®

This feature pauses chair movement when disrupting impact is detected.

6. PELVIC TILT

Tilting of the seat section assists with proper exam positioning.

7. FRONT STORAGE DRAWER

The additional storage gives you room to store supplies at the point of care.

8. BARRIER-FREE LOW CHAIR HEIGHT

The 46 cm low seat height makes accessing the exam chair easier for all patients and complies with US ADA recommendations



Midmark 627 shown with
Patient Support Rails+ accessory.

*Features available on the 627 model only.